

# Chromatographic Figures of Merit



**Agilent Technologies**


## In This Section, We Will Discuss:

- The relevance and use of capacity factor
- The calculation of efficiency or plate number
- How to calculate the resolution between two chromatographic peaks
- Peak symmetry

**NOTE: You will need a hard copy of the chromatogram and report from the last laboratory exercise.**

# Capacity Factor and Relative Retention

- **Capacity factor** is characteristic of a specific compound at a given mobile phase composition, temperature, and column type.
- **Capacity factor** is equal to the number of moles in the stationary phase divided by the number of moles in the mobile phase.



Capacity Factor

$$k' = \frac{t_R - t_o}{t_o}$$

Relative Retention

$$t'_R = t_R - t_o$$



# Calculate Capacity Factor




Capacity Factor is independent of flow rate  
making day-to-day fluctuations less troublesome



Calculate the capacity factor for the  
third major peak in your chromatogram.




# Resolution



For equal peak areas, R of 1.5  
gives baseline separation



# Calculate Resolution



R - resolution

$t_{RB}$  - retention time of component B

$t_{RA}$  - retention time of component A

w - width at base of peak

$w_{1/2}$  - width at half-height

$$R = 2 \left( \frac{t_{RB} - t_{RA}}{W_A + W_B} \right) \quad R = 1.176 \left( \frac{t_{RB} - t_{RA}}{W_{1/2A} + W_{1/2B}} \right)$$




Calculate the resolution between the 2nd and 3rd chromatographic peaks.




Agilent Technologies

# Efficiency



# Calculate Efficiency



Calculate the efficiency of the fourth peak.

$$N = 16 \left( \frac{t_R}{W_B} \right)^2 = 5.54 \left( \frac{t_R}{W_{1/2}} \right)^2 = 2 \Pi \left( \frac{h_{ptr}}{A} \right)^2$$

$$HETP = \frac{L}{N}$$

*N: Efficiency*

*HETP: Height Equivalent to a Theoretical Plate*

*L: Column Length*



*hp: Peak Height*

*A: Peak Area*





# Peak Symmetry



# Chromatographic Figures of Merit

These calculations can be done automatically by your ChemStation:

- Select a System Suitability (Performance) report style
- Perform a Sequence Summary Report

